δ^+ PolyK δ^-

PK-C30kV-100C

Corona Poling System for Research

- Poling Voltage: up to -30 kV
- Temperature: up to 100 ° C
- Size: 100 mm x 100 mm

PolyK Technologies State College, PA, USA

www.piezopvdf.com

sale@polyk-lab.com

Corona Poling System © PolyK Technologies

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
-30 kV Output Cable

Please do NOT touch the connectors during poling!

The grounding connectors can be wrapped with electrical insulation tape.

Grounding Cable
from poling
station to building
electrical
grounding

HV Cable from
Power Supply
(female) and
from poling
station (male)

The grounding connectors
can be wrapped with
electrical insulation tape.

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

www.piezopvdf.com

sale@polyk-lab.com

Thank you for choosing our ferroelectric test instruments

- Specialized in high voltage polymer dielectric, piezoelectric, pyroelectric, and smart materials, test instrument, and applications
- Provide various PVDF (different Mw) homopolymer resin, P(VDF-TrFE) (6 compositions), P(VDF-TrFE-CTFE), P(VDF-TrFE-CFE), P(VDF-HFP), P(VDF-CTFE), P(VDF-TFE), P(TFE-HFP-VDF), etc.
- Thin dielectric film, piezoelectric and pyroelectric PVDF and PVDF-TrFE film, uniaxially oriented or biaxially oriented.
- Ultrathin high temperature polymer dielectric film, PET, BOPP, PEN, PPS, PC, 0.9 μm to 10 μm
- Special low-cost high voltage dielectric and ferroelectric test systems.

For help, please contact sale@polyk-lab.com

Corona Poling System © PolyK Technologies www.piezopvdf.com

13